

## Hospital ICU Bed



**Product Categories:** [Hospital Equipment & Furniture](#)

**Product Tags:** [hospital](#), [medical](#)

**Product Page:** <https://www.labappara.com/product/hospital-icu-bed/>

### Product Description

Hospital ICU Bed - Hospital Furniture & Equipment

Top of CRCA sheet in four sections perforated uniformly. Trolley type base & Main work precise steel tubes. Base Mounted on 125 Mm dia swiveling castors two with Brakes.

The backrest and knee rest also Trendelenburg and reverse movement controlled by the actuator.

The level height of bed can be facilitating to close observation by the medical staff. Stainless Steel bows with laminated panels easily removable for access to the patient.

Fitting IV Rod & and protective side railing.

**Overall Size** - 208L X 90W X 60/98H cm