

Hospital Semi Flower Bed



Product Categories: [Hospital Equipment & Furniture](#)

Product Tags: [hospital](#), [medical](#)

Product Page: <https://www.labappara.com/product/hospital-semi-flower-bed/>

Product Description

Hospital Semi fowler Bed - Hospital Furniture

Mechanically operated two sectional, uniformly perforated CRCA sheet top beds. Firm commitment to quality, we are presenting for our clients a vast range of Semi Fowler Bed. The offered fowler bed is precisely designed at our production unit by utilizing the superior quality components in compliance with defined guidelines. This fowler bed has single screw mechanism which helps in adjusting the bed for backrest.

Feature -

Strong "H" Shaped legs on rubber stumps.

One Fix lever handles at foot end to obtain head raise upto 80deg.

Overall Size - 208L X 92W X 60H