

Human molar Teeth



Product Categories: [Biology Lab](#), [Educational equipments](#), [Educational Models](#), [Skeleton Models](#)

Product Tags: [Educational Model](#)

Product Page: <https://www.labappara.com/product/human-molar-teeth/>

Product Description

Human Molar Teeth- Scientific Educational Model

Labappara Human Molar Teeth model illustrated the structure of the molar teeth , also show the capillary vessel and nerve. humans, the molar teeth have either four or five cusps. Human have twelve molars, in four groups of three at the back of the mouth. The third, molar in each group is called a wisdom tooth.

This Model Depicts An Upper Triple- Root Molar and Separates Into 6 Parts.

It Features a Longitudinal Section Through the Crown, Two Roots and a Pulp Cavity.

Contains Removable Pulp and Three Tooth Inserts with Different Stages of Advanced Caries. On Stand

Specifications:

Height 11cm.

Width 15cm.

Thickness 4cm.

N.W. : 0.9 Kg.

Product data have been exported from - [Labappara scientific instrument](#)

Export date: Wed Nov 27 22:38:21 2024 / +0000 GMT