

Huygens Eyepiece 10X, 15X



Product Categories: [Eye Piece](#), [Microscope](#), [Microscope Assecories](#), [Optical Instruments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/huygens-eye-piece-10x-15x/>

Product Description

Huygens Eyepiece 10X, 15X - Scientific Equipment

Huygens design consists of two Plano-convex lenses positioned with both convex faces toward the object.

Designed for use with low power achromatic objectives.

Our entire line of eyepieces will fit any international standard microscope with a 23 mm inner diameter eyepiece tube.