

# Huygens Eyepiece 5X



**Product Categories:** [Eye Piece](#), [Microscope](#), [Microscope Assecories](#), [Optical Instruments](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/huygens-eye-piece-5x/>

## Product Description

Huygens Eyepiece 5X - Scientific Equipment

Huygens design consists of two Plano-convex lenses positioned with both convex faces toward the object.

Designed for use with low power achromatic objectives.

Our entire line of eyepieces will fit any international standard microscope with a 23 mm inner diameter eyepiece tube.