

Hydrostatic Balance (Brass)



Product Categories: [Balance & Scale](#), [Chemistry Equipments](#)

Product Tags: [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/hydrostatic-balance-brass/>

Product Description

Hydrostatic Balance (Brass) - Laboratory Equipment

These are general purpose balances for students.

The hangers have double hooks for specific gravity experiments, and the pans are detachable.

Feature -

Finished in durable black enamel with chrome plated.

Fitting and mounted on nicely polished heavy wooden base.

Provided with leveling screws.

The balance is mounted on a rectangular wooden polished base.