

HYPSONETER COPPER



Product Categories: [Heat And Thermodynamics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/hypsoneter-copper/>

Product Description

Hypsoneter Copper - Physics Equipment

For determining the upper fixed point of thermometers and boiling points at Barometer pressure.

Comprising a Copper boiling vessel with a steam jacketed in a tube.

A lid with a tubule is provided, which holds the thermometer under test in a cork

A side outlet draining tube and a glass Mano metric tube mounted in tubular with rubber bung.

Supplied without thermometer.

SIZE:- 8", 10", 12"