

## Immersion Thermometer



**Product Categories:** [Chemistry Equipments](#), [Thermometer & Hygrometer](#)

**Product Tags:** [chemistry lab instruments](#), [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/immersion-thermometer/>

### Product Description

Immersion Thermometer - Laboratory Equipment

Thermometers manufactured for total immersion indicate an accurate temperature reading

When the thermometer is immersed to the level of the liquid in the capillary column.

To begin, immerse the thermometer to the top of the liquid column as the temperature rises or falls adjust the thermometer until an accurate reading is achieved.

**Available in sizes:**

1100 C X 1/1 24", 36"

2500 C X 1/1 24", 36"

3600 C X 1/1 24", 36"