

Inclined Medical Microscope(Deluxe Model)



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page:

<https://www.labappara.com/product/inclined-medical-microscope-delux-model/>

Product Description

Inclined Medical Microscope(Deluxe Model) - Scientific Equipment

Deluxe Model heavy body fitted 6 V-20 W halogen lamp workable on 220 Volt AC provided WITH Coaxial mechanical stage in thick plywood box. Optical combination. A large stage fitted with COAXIAL MECHANICAL STAGE having low drive co- axial knobs for X-Y movement of slide.

A built in base 6 V-20 W halogen light illumination with Light Intensity regulator. Provision for easy replacement of Halogen bulb.

Superior Quadruple revolving Nosepiece with positive click stops.

Observation head is 45° MONOCULAR INCLINED TUBE rotatable to 360°.

Low position Coarse & calibrated Fine Focusing controls.

Optical Combination :-

Objectives : Achromatic 5 x, 10 x, 40 x/45 x (SL) & 100 x SL Oil Immersion (SL) (Self Parafoalised)

Eye Pieces: Coated WF 10 x & Huy 5 x.

Packed in a Styrofoam box inside an attractive Sun-mica clad plywood box.

Features:

Crystal clear images

Trouble-free functioning

Compact design