

Inclined Medical Microscope(Economy Model)



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page:

<https://www.labappara.com/product/inclined-medical-microscope-economy-model/>

Product Description

Inclined Medical Microscope(Economy Model) - Laboratory Instrument
Economy Model heavy body fitted 6 V-20 W halogen lamp workable on 220 Volt AC provided WITH Co-axial mechanical stage in thick plywood box. Optical combination.

All aluminum pressure die casted powder coated body gives it a nice look, shape and is available in various desired shades. Inclined Microscope has the quadruple ball bearing nose-piece with rotate able observation head 360° and 45° inclined

head with all anti fungus coated prisms. The separate knobs provided for cores and fine motion have best count of 0.002 mm to focus system. The various set of Eyepiece & Objectives can be fitted like 4 x, 10 x, 40 x, and 100 x and paired set of 10 x wide field eyepieces. It is illuminated by 6 V/20 W halogen lamp needs 220 V & 50 Hz supply.

Features:

Crystal clear images

Trouble-free functioning

Compact design

Eye piece 5x, 10x (Hug)

Objective 5x, 10x, 45x (SL), 100x Oil Immersion (SL) (Self Parafofocalised)