

Inclined Monocular Microscope



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page:

<https://www.labappara.com/product/inclined-monocular-microscope/>

Product Description

Inclined Monocular Microscope - Scientific Equipment

Inclined microscope with inclined tube rotatable to 360° disc diaphragm and built in condenser. Separate coarse motion & highly sensitive slow motion working on ball bearing guide-ways. Quadruple revolving nose-piece, attachable mechanical stage. Illumination by 6 V-20 W halogen bulb with solid state transformer. Separate 55 mm dia reflector for day light illumination. Complete with lens cleaning paper,

immersion oil & other standard accessories. Supplied in an attractive & thick plywood box.

Features:

High magnifying power

Provides accurate results

Easy to use

Achromatic objective: 4 x/5 x, 10 x, 40 x/45 x & 100 x Oil immersion. Eyepiece: 10 x WE or 15 x WE

Magnification: 40 x to 1500 x.