

## Incubator Universal



**Product Categories:** [Incubators](#), [Lab Heating- Equipment](#), [Lab Incubators](#), [Laboratory Equipments](#)

**Product Tags:** [biology lab equipment](#)

**Product Page:** <https://www.labappara.com/product/incubator-universal/>

### Product Description

Incubator Universal - Biology Instrument

In biology an incubator universal is a device used to grow and maintain microbiology culture the carbon dioxide (CO<sub>2</sub>) and oxygen content of the atmosphere inside.

**Description:**

Temperature up to 250°C. heating elements are placed in ribs at the bottom and sides.

Double walled inside Stainless Steel.

Outside Mild Steel painted in epoxy powder coating/stainless steel with perforated

adjustable shelves, two indicators and on/off switch.

To work on 220/230 volts A.C. with OR without one air circulating fan & digital temperature controller cum indicator.

**Features:**

High functionality

Energy efficient

Minimum maintenance

**Aluminum chamber-**

300 x 300 x 300 MM (12" x 12" x 12")

350 x 350 x 350 MM (14" x 14" x 14")

450 x 450 x 450 MM (18" x 18" x 18")

600 x 600 x 600 MM (24" x 24" x 24")

**Stainless Steel chamber-**

300 x 300 x 300 MM (12" x 12" x 12")

350 x 350 x 350 MM (14" x 14" x 14")

450 x 450 x 450 MM (18" x 18" x 18")

600 x 600 x 600 MM (24" x 24" x 24")