

Incubators Stainless Steel (B.O.D.)



Product Categories: [Incubators](#), [Lab Heating- Equipment](#), [Lab Incubators](#), [Laboratory Equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page:

<https://www.labappara.com/product/incubators-b-o-d-stainless-steel-2/>

Product Description

INCUBATORS (B.O.D.) STAINLESS STEEL - Laboratory Equipment

Temp range 50 C to 500 C with accuracy + 1 C.

Double walled, inside anodized aluminum Stainless Steel and outside M.S.

Power Coated to work on 220/230 volt A.C. 50 h

The low temperature incubator is made of double walled sheet finished with epoxy based powder coating.

Inner chamber is made of highly polished stainless steel / anodized aluminum.

Two doors are provided of INCUBATORS (B.O.D.) STAINLESS STEEL.

Inner door is made of transparent acrylic for inspecting specimens.

Outer door is insulated and is fitted with magnetic tape with lock and key.

Provided with two air circulating fans, high performance compressor , cooling coils, heating elements, digital Temperature control and display is also fitted with a door operated illumination and a castor wheel for easy mobility.

Inner Chamber Size In mm (W X D X H) Capacity In Cu.Ft

Capacity

455 x 410 x 610

4

112 ltrs

505 x 415 x 830

6.1

171 ltrs

565 x 565 x 865

10

280ltrs

650 x 580 x 900

12

336 ltrs

700 x 650 x 900

15

420 ltrs

Optional:

Time 0-2 Hours

Digital Temp, Indicator-cum-controller with timer.

Microprocessor PID digital temp, Indicator-cum-controller with timer.

Automatic voltage Stabilizer 3-5KVA (recommended must).

Co2 air mixer nozzle is Provided