

Infant Nursing Simulator



Product Categories: [Biology Lab](#), [Educational equipments](#), [Paediatrics](#), [Gynae & Obstetric Model](#), [Reproductive Models](#)

Product Tags: [biology lab](#), [fibre model](#), [gynae model](#), [paediatrics](#), [science model](#)

Product Page: <https://www.labappara.com/product/infant-nursing-simulator/>

Product Description

Infant Nursing Simulator- Scientific Educational Model Specification

Realistic fontanel, coronal suture and parietal suture.

Place nasal and oral gastric tubes.

Temperature measurement.

Umbilicus nursing.

Realistic testis in the scrotum.

Training of bathing, lactation, lavage, etc.

Features :

1. Head, neck and limbs can be moved freely. The chest skin can be replaced;
2. Cephalic venipuncture;
3. Umbilical venipuncture;
4. Oral and nasal intubation nursing;
5. Care of gastrolavage, enema and urethral catheterization;

6. Ostomy;
7. Deltoids and buttocks intramuscular injection;
8. Holistic nursing; sponge bath, replacing clother, lactation and diaper changing.