

Bell Jars Round Bottom



Product Categories: [Desiccators Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [glassware](#)

Product Page: <https://www.labappara.com/product/jars-bell-round-bottom/>

Product Description

Bell Jars, Round Bottom - Laboratory Glassware Equipment

A bell jar is a piece of laboratory equipment used to contain vacuum. It is similar in shape to a bell.

These are made of thin glass, with more care being taken regarding their optical clarity .

These jars come with properties like transparent & superior heat-resistance.

Available in O.D. x Height (mm. approx.) - 200 x 350 , 300 x 400, 300 x 450,