


## Joule's Electrical Calorimeter



**Product Categories:** [Analytical Lab Instruments](#), [Calorimeters](#)

**Product Tags:** [science lab equipment](#)

**Product Page:**

<https://www.labappara.com/product/joules-electrical-calorimeter-seasoned-wood/>

### Product Description

Joule's Electrical Colorimeter - Analytical Laboratory Equipment

The unit consists of a Copper calorimeter with outer vessel of well-polished wood. Fitted with Bakelite top with hole for thermometer and stirrer. Two terminals are connected to a coil of wire for passing current.

#### **Features:**

Ideal for using college or research laboratory

Easy to read precisely spaced markings

High sensitivity and accurate temperature indication

Compact design

Light-weight

Low maintenance

Complete with stirrer.

**Available in size - 3"x 2" , 4"x 3"**