


Joule's Electrical Colorimeter (With Teak Wood)



Product Categories: [Analytical Lab Instruments](#), [Calorimeters](#)

Product Tags: [science lab equipment](#)

Product Page: <https://www.labappara.com/product/joules-electrical-calorimeter-teak-wood/>

Product Description

Joule's Electrical Colorimeter (With Teak Wood) - Analytical Equipment

Consists of a copper vessel With Teak Box of wood fitted with Ebonics top with holder for thermometer.

The terminal are connected to a coil of wire for 3 amps.

Complete with stirrer.

Available in size - 3"x 2" , 4"x 3"