

Kelvin Bridge With Head Phane & Audio Oscillator 1 KHz



Product Categories: [Electrical Bridge's](#), [Electrical Bridge's](#), [Electronics](#), [Engineering Equipment](#)

Product Page:

<https://www.labappara.com/product/kelvin-bridge-with-head-phane-audio-oscillator-1-khz/>

Product Description

Kelvin Bridge With Head Phane & Audio Oscillator 1 KHz

Features:

User-friendly operation

Longer serving life

Compact design

Other Details:

Dedicated trainer board for each bridge.

Completely self-contained stand-alone unit.

Demonstrates the principle and working of ac bridges for measurement of resistance.

Supply required 230v, 50 hz ac.

Built-in ic based dc regulated power supply with short circuit protection and led indication for supply "on".

Built-in 1khz oscillator.

Built-in imbalance amplifier.

Head-Phone set for sensitive detection.

Built-in variable arms and multiplier.

Multi-coloured test points at various stages in the circuit to observe waveforms and voltages.

Set of patch cords.

Housed in an elegant cabinet with a well spread intelligently designed circuit layout on the front panel.

Strongly supported by a comprehensive manual complete with theory and operating details.