

Keyboard Interface



Product Categories: [Electronics](#), [Engineering Equipment](#), [Interfacing Cards](#), [Interfacing Cards](#), [Microprocessor & Micro-controller Lab](#), [Microprocessor & Micro-controller Lab](#)

Product Page: <https://www.labappara.com/product/keyboard-interface/>

Product Description

Keyboard Interface

A calculator style kb consisting of the keys 0 to 9, *, +, -, /, c ce and two user definable keys, allowing the user to study a number of techniques used in key board interfacing like software debouncing, two key lockout and keyboard encoding and pausing.