

# Clinostat



**Product Categories:** [Plant Physiology Apparatus](#)

**Product Page:** <https://www.labappara.com/product/klinostat/>

## Product Description

Clinostat - Plant Physiology Equipment

Fitted on a heavy metal stand. It can run 30 hrs. in vertical or horizontal position and to carry weights with clock movement.

Clinostat with new clock movement having start and stop system.

And then fitted on heavy cast iron base.

Clinostat a plant carrying cup is also provided.

A Clinostat is mainly a device that makes use of rotations for negating the effects of gravitational pull on plant growth and development.

Clinostat is used to carry weights based on heavy metal stands in a clockwise movement.

It can also run continuously for 30 hours in vertical or horizontal position.