

## KUNDT'S TUBE



**Product Categories:** [Heat And Thermodynamics](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/kundts-tube/>

### Product Description

Kundt's Tube - Physics Equipment

Used To determine the standing longitudinal waves in column of air and velocity of sound in air.

It comprises of borosilicate glass tube approx. 1.2 meter Length and 34 mm diameter.

There is a tightly fitting but adjustable piston on a 390 mm long rod at one end of the glass tube.

A loose fitting piston attached to a 1.2 meter long sounding rod at the other end of the glass tube.

Fitted on a wooden base board.