

Kymograph



Product Categories: [Pharmaceutical Equipments](#)

Product Tags: [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/kymograph-2/>

Product Description

Kymograph - Pharmaceutical Equipment

Cylinder 6" x 11" with crown wheel engaging shaft. A spring loaded device allows turning by hand if necessary.

Moving in slot in the bed allowing variation in paper length for 60" to 66".

Rigid cast iron stand works up & down by a large screw.

Fine adjustment ensuring return of outer arm to exactly the same position after change of another paper.

Kymograph (which means wave writer) is a device that gives a graphical

representation of spatial position over time in which a spatial axis represents time. It basically consists of a revolving drum wrapped with a sheet of paper on which a moves back and forth recording perceived changes of phenomena such as motion or pressure.

Features:

Precisely engineered

Sturdiness

Easy to adjust

Accessories supplied with equipment (optional) :

Double Time Maker.

Mercury Manometer.

Separate smoking & Varnishing Apparatus.

Smoking Burner.