


Kymograph



Product Categories: [Physiology Apparatus](#)

Product Page: <https://www.labappara.com/product/kymograph/>

Product Description

Kymograph - Laboratory Equipment

Electric recording drum 4/8 speed gear system or 8 speed direct.

Electrically driven drum suitable for general experiment working.

Kymograph with S.S spindle on ball bearing moment, aluminum cylinder 6" dia, robust gear clutch.

The cylinder is raised or lowered and is locked in position by a small lever.

Accessories for Kymograph -

Accessories for kymograph set of 7 for conducting simple muscle twitch experiment viz.

Accessories for adjustable stand induction coil, simple electrode.

Tuning fork N-100, weight set & power supply.

Kymograph is in glazed paper sheet.

Features:

Requires less maintenance

Robust design

Easy to operate

Extensively used

Durable

Highly reliable