

Laboratory Spatula



Product Categories: [Connectors And Stirrers](#), [Laboratory Plastic Ware](#)

Product Tags: [Laboratory Equipments](#), [Laboratory plastic ware](#)

Product Page: <https://www.labappara.com/product/laboratory-spatula/>

Product Description

Laboratory Spatula - Scientific Laboratory Equipment

Moulded in Polypropylene, Spatula consists of a spoon on one side while the other end is a flat construction with a knife edge. This edge is helpful in breaking up crystals of different salts while the spoon is used to draw salts from the container

Particulars	Packing
150 mm	12 Pcs
200 mm	12 Pcs