

Labour Table



Product Categories: [Hospital Equipment & Furniture](#)

Product Tags: [hospital](#), [medical](#)

Product Page: <https://www.labappara.com/product/labour-table/>

Product Description

Labour Table - Hospital Equipment & Furniture

Constructed from MS tubular steel material. Adjustable head and leg section. Top is upholstered and covered with water proof synthetic material.

Complete with SS bowl 32cm dia.

Standard padded leg holder.

Finish in epoxy power coated.

Overall Size - 1830 L X 550 W X 800 Hmm