

LCR Impedance Apparatus



Product Categories: [Characteristics And Application](#), [Characteristics And Application](#), [Electronics](#), [Engineering Equipment](#)

Product Page: <https://www.labappara.com/product/lcr-impedance-apparatus/>

Product Description

LCR Impedance Apparatus

Features :

Instrument comprises of variable AC Power Supply 0-15 VAC, 2 AC moving coil meters to measure voltage & current, 2 Resistance, 2 Capacitors & 1 Inductance connected inside & connections brought out at terminals