

## Litmus Paper



**Product Categories:** [Chemistry Equipments](#), [General Chemistry Labware](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [science lab equipment](#)

**Product Page:** <https://www.labappara.com/product/litmus-paper/>

### Product Description

Litmus Paper - Chemistry Laboratory Equipment

Paper stained with litmus which is used to indicate the acidity or alkalinity of a substance.

Blue litmus paper turns red under acid condition.

Red litmus paper turns blue under alkaline conditions

Packing Unit: 100 LVS(5 books X 20LVS) 200 LVS(10 books X 20LVS)

Storage: in Acid & Alkali free zon

Stability: 5 Years

Details: Made from imported filter paper

pH indicator paper turns red under acidic conditions. pH range 4.3 - 6.8. 100 strips.

Blue litmus paper Blue , pH range 8-5 . Changed color from blue to red in acid.