

## Logic Controller Interface



**Product Categories:** [Electronics](#), [Engineering Equipment](#), [Interfacing Cards](#), [Interfacing Cards](#), [Microprocessor & Micro-controller Lab](#), [Microprocessor & Micro-controller Lab](#)

**Product Page:** <https://www.labappara.com/product/logic-controller-interface/>

### Product Description

#### Logic Controller Interface

Logic Controller Interface is worldwide used for making multiple arrangements in digital as well as analog inputs and outputs for extending temperature range.

Logic Controller Interface Card The logic controller interface provides the user with eight TTL buffered outputs and eight TTL buffered inputs. The logic state of each input and output is indicated by LED. This module can be effectively used to teach 8255 input & output modes . It can also simulate Ladder networks..

#### Features:

Excellent functionality

Shock proof

Superior performance