

Losses In Pipe Fitting Apparatus



Product Categories: [Engineering Equipment](#)

Product Tags: [Engineering Equipment](#)

Product Page: <https://www.labappara.com/product/losses-pipe-fitting-apparatus/>

Product Description

Losses In Pipe Fitting Apparatus- Engineering Equipment

The unit provide study of losses of pipe in various pipe fitting. A unit consists of Major and minor losses associated with pipe flow in piping networks helps in determining the pumping power requirements, material and fittings selection. A good understanding of such losses helps engineers in designing optimum fluid distribution systems, process plants etc. This apparatus is designed to introduce students to minor flow losses (Fittings) in Pipes

Feature

Block Type Acrylic Differential Manometer (250-0-250 mm)

Number of fittings : 4 Nos

Flow control valve to change Discharge.

Transparent pizometer tubes with scale to measure measuring tank discharge.

Measuring tank Dimensions- 295 X 345 X 345 mm

Sump tank Dimensions- 900 X 350 X 350 mm