

## Magnesium Metal Ribbon



**Product Categories:** [Chemistry Equipments](#), [General Chemistry Labware](#)

**Product Tags:** [chemistry lab instruments](#), [science lab equipment](#)

**Product Page:** <https://www.labappara.com/product/magnesium-metal-ribbon/>

### Product Description

Magnesium Metal Ribbon - Chemistry Laboratory Equipment

A light, silvery-white, moderately hard metallic element that in **ribbon** or powder form burns with a brilliant white flame.

Features:

Highly effective

Easily react with acids

Reliable result

Specifications:

Formula: Mg

M.W: M.W. 24.31

Min.assay: 99.50%

Fineness 99.95%

Weight 25 Gm

Length Approx 21 m (73 ft)

Width 4 mm

Thick 0.2 mm