

# Magnetic Amplifier



**Product Categories:** [Control Lab](#), [Control Lab](#), [Electronics](#), [Engineering Equipment](#)  
**Product Page:** <https://www.labappara.com/product/magnetic-amplifier/>

## Product Description

Magnetic Amplifier

### Technical Specifications

Magnetic Amplifier : High gain, vacuum impregnated with twin storable reactors.

Reactor : Storable, made of torroidal core.

Core : Torroidal, cold rolled grain oriental core with low loss & high permeability.

Windings : 1. Load (gate) windings wound dividually on each saturable eactors.

2. Control winding.

3. Bias winding.

4. Feedback winding, control bias and feedback windings are wound on the two saturable reactors combined together.

Rectifier : Set of diodes.

Transformer : Power transformer

Load Capacity : 100W ( max ) facility for connecting external load.