

Magnetic Field Apparatus



Product Categories: [Educational equipments](#), [Electromagnetism](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/magnetic-field-apparatus/>

Product Description

Magnetic Field Apparatus - Physics Equipment

Description :

Magnetic Field Apparatus is used to demonstrate the production of magnetic field by current flowing in a solenoid. A small wooden table supports a Bakelite ring with three coils of enameled copper wire of 5, 10, and 15 turns respectively connected to the terminals at the edge of the base.

Features :

Sensitive

Compact

Flexible

Adjustable coil spacing

Easy to use

Stable

Reliable