

## Magnetic Needle Stand (Plastic)



**Product Categories:** [Educational equipments](#), [Electromagnetism](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:**

<https://www.labappara.com/product/magnetic-needles-stand-plastic/>

### Product Description

Magnetic Needle Stand (Plastic) - Physics Equipment

Magnetic Needle Stand is made up of high quality plastic and is used to mount magnetic needle.

Magnetic Needle mounted on Stand is built to indicate the earth's magnetic fields or position of North and South Pole approximately.