

Magnetic Needles



Product Categories: [Educational equipments](#), [Electromagnetism](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/magnetic-needles/>

Product Description

Magnetic Needles - Physics Equipment

Made up of Special Steel, Highly Magnetized With brass bearing for pivoting.

Magnetic Needle Jewel Bearing of Magnetic Needle Jewel Bearing from India. We bulk supply Magnetic Needle Jewel Bearing for tenders.

Available in size -50 mm, 62 mm, 75 mm, 100 mm.