

## Magnifiers (Bakelite Handle)



**Product Categories:** [Educational equipments](#), [Light & Optics](#), [Magnifier](#), [Optical Instruments](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/magnifiers-bakelite-handle-3/>

### Product Description

#### Magnifiers (Bakelite Handle) - Physics Optical Instrument

These hand held Bakelite Handle Magnifier having Spherical Double Convex glass lens mounted in chromium plated brass ring with metal knob and Bakelite handle looks very attractive and specially manufactured for export purpose.

The hand-held magnifier design involves a handle that is attached to a lens through which the user will view an object - whether it is text or an image Hand-held magnifiers are the one of the simplest type of magnifiers to use and are also the most familiar. Nearly everyone has used a hand-held magnifier at some point in time.

<b>Diameter</b>	<b>Focal length</b>
-----------------	---------------------

75mm	20-100cm
------	----------

100mm	20-100cm
-------	----------