

## Magnifiers (C.P Handle brass C.P Ring)



**Product Categories:** [Educational equipments](#), [Light & Optics](#), [Magnifier](#), [Optical Instruments](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:**

<https://www.labappara.com/product/magnifiers-c-p-handle-brass-c-p-ring/>

### Product Description

Magnifiers (C.P Handle brass C.P Ring) - Physics Optical Equipment

A magnifier is a device used for magnification. A magnifying glass (called a hand lens in laboratory contexts) is a convex lens that is used to produce a magnified image of an object. The lens is mounted in a frame with a handle . A magnifying glass can also be used to focus light, such as to concentrate the sun's radiation to create a hot spot at the focus for fire starting.

Magnifier can also refer to: Magnifying glass, an optical device for magnification  
Screen magnifier, software.

<b>Diameter</b>	<b>Focal length</b>
-----------------	---------------------

37mm	20-100cm
------	----------

50mm	20-100cm
------	----------

75mm	20-100cm
100mm	20-100cm
150mm	20-100cm