

Mass Transfer With & Without Chemical Reaction (Solid -Liquid)



Product Categories: [Engineering Equipment](#), [Mass Transfer Lab](#)

Product Page: <https://www.labappara.com/product/mass-transfer-without-chemical-reaction-solid-liquid/>

Product Description

Mass Transfer With & Without Chemical Reaction (Solid -Liquid)

Description:

The setup consists of a Reactor, in which solid cylinder of Benzoic Acid is placed and is stirred by means of a variable speed mixer. Baffles are also provided. The Reactor is made of Stainless Steel. Peltizer is provided to make the solid cylinder of Benzoic Acid. Constant temperature water bath arrangement is provided to conduct the experiment at various temperatures.