

Maternal & Neonatal Delivery Simulator



Product Categories: [Biology Lab](#), [Educational equipments](#), [Paediatrics](#), [Gynae & Obstetric Model](#), [Reproductive Models](#)

Product Tags: [biology lab](#), [fibre model](#), [gynae model](#), [paediatrics](#), [science model](#)

Product Page:

<https://www.labappara.com/product/maternal-neonatal-delivery-simulator/>

Product Description

Maternal & Neonatal Delivery Emergency Simulator: Scientific

Educational Model Maternal manikin features:

The mechanical gear device is equipped with two mechanical adapters for connecting stimulative delivery fetus; fetus and adapter, the adapter and adapter, adapter and gear device are all fixed with elastic clamp device.

The delivery process and fetus heart sound controller can suspend, initialize, start and continue

The delivery process; delivery speed is according to requirement, have 1 to 4 steps;

Fetal heart sound auscultation: can set fetal heart sound frequency and volume, heart rate is adjustable between 70-180;

Can simulator head birth, delivery, narrow birth canal, cord around neck, fore-lying placenta, etc;

Equipped with lifting "Soft pad" for practicing Leopold maneuver;

Can measure fetal head descent and cervical dilation

Can simulate multiple placenta position

Can set venous channel in maternal arm for giving administration and nutrition

Vulvae suture practice modules, including left lower, medial and right lower three incision position;

Tracheal incubation training

Manually simulate carotid artery pulse.

Neonatal manikin features:

Phosphoric function

Nursing features: eye irrigation and drops, neonatal washing and bandaging;

Oral and nasal incubation, infant sputum suction, tracheal incubation and gastric lavage;

Infant umbilical cord nursing, scalp vein puncture, arm venous puncture, there is obvious feeling for venous puncture, and blood flash back indicated correct puncture

Components:

Advance Delivery and Maternal and Neonatal Emergency Simulator System.

Neonatal simulator system (for deliver use)

Neonatal Emergency within Nursing model

Automated delivery and fetal heart sound controller

Adult CPR electronic display

BP Measurement Trainer

Elasticity Cervix

Cervix examination modules 5 pcs/set

Saturate perinealis model 3 pcs/set

48 hours postpartum uterus

Wooden Stethoscope,use for fetus heart

Normal placenta (1 set)

Neonatal airway cannual

Blood pressure Meter (with cuff)

Product data have been exported from - [Labappara scientific instrument](#)

Export date: Wed Nov 27 16:45:07 2024 / +0000 GMT