

## Maximum & Minimum Thermometer (Duplex Box)



**Product Categories:** [Chemistry Equipments](#), [Thermometer & Hygrometer](#)

**Product Tags:** [chemistry lab instruments](#), [physics lab instruments](#)

**Product Page:**

<https://www.labappara.com/product/maximum-minimum-thermometer-duplex-box/>

### Product Description

Maximum & Minimum Thermometer (Duplex Box) - Laboratory Equipment

Comprising two -10 to - 50 °C x 1 °C thermometers fitted on scale. One thermometer bulb is covered with muslin which kept moistened being connected to a water reservoir. Mounted in a case with hinged cover.

### Features:

Impeccable finish

High accuracy

Precisely designed