


MCT Twin Rack



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/mct-twin-rack/>

Product Description

MCT Twin Rack - Scientific Laboratory Plastic ware

Labappara provides this Twin Purpose Rack for both micro centrifuge tubes. Each side has holes for MCT's respectively. The transparent top provides an excellent view as well as a protection cover for the tubes stored. Each face for the base has a unique alpha-numeric index for easy identification of the tube. Top & base of this rack are made of tough & transparent polycarbonate while the middle portion is made of polypropylene. The tube holding spaces are uniformly placed so as to facilitate the use of multi-pipettes. Sturdy polypropylene rack accommodates individual 0.2 ml tubes.