


# Measuring Cylinder (Hexagonal Base)



**Product Categories:** [Laboratory Plastic Ware](#), [Liquid-Measurement](#)

**Product Tags:** [Laboratory Equipments](#), [Laboratory plastic ware](#)

**Product Page:**

<https://www.labappara.com/product/measuring-cylinder-hexagonal-base/>

## Product Description

Measuring Cylinder (Hexagonal Base) - Scientific Laboratory Equipment  
Base measuring cylinders are now available with hexagonal base also. The hexagonal base provides excellent stability to the cylinder where as the tapered spout works really well while pouring any agent.

Particulars	Least Count	Packing
10ml	0.2 ml	24 Pcs
25ml	0.5 ml	12 Pcs
50ml	1 ml	12 Pcs
100ml	1 ml	12 Pcs
250ml	2 ml	6 Pcs
500ml	5 ml	4 Pcs
1000ml	10ml	3 Pcs