

## Measuring Cylinders (Round Base)



**Product Categories:** [Laboratory Glassware](#), [Pipettes](#), [Volumetric Flasks And Measuring Cylinder](#)

**Product Tags:** [glassware](#), [Laboratory Equipment](#)

**Product Page:**

<https://www.labappara.com/product/measuring-cylinders-round-base/>

### Product Description

Measuring Cylinders (Round Base) - Laboratory Glassware Equipment

Measuring Cylinders are made from borosilicate glass. A graduated cylinder, measuring cylinder or mixing cylinder is a common piece of laboratory equipment used to measure the volume of a liquid. It has a narrow cylindrical shape. Each marked line on the graduated cylinder represents the amount of liquid that has been measured.

#### **Applications:**

Used to measure the volume of a liquid

#### **Features:**

Temperature resistance

Crack resistance

Precisely designed

Capacity - 5 ml, 10 ml, 25 ml, 50 ml, 100 ml, 250 ml, 500 ml, 1000 ml, 2000 ml