


Measuring Cylinders



Product Categories: [Laboratory Glassware](#), [Pipettes](#) , [Volumetric Flasks And Measuring Cylinder](#)

Product Tags: [glassware](#)

Product Page: <https://www.labappara.com/product/measuring-cylinders/>

Product Description

Measuring Cylinders - Laboratory Glassware Equipment

Measuring Cylinders are made from borosilicate glass. A graduated cylinder, measuring cylinder or mixing cylinder is a common piece of laboratory equipment used to measure the volume of a liquid. It has a narrow cylindrical shape. Each marked line on the graduated cylinder represents the amount of liquid that has been measured.

Applications:

Used to measure the volume of a liquid

Features:

Temperature resistance

Crack resistance

Precisely designed

Capacity - 5 ml, 10 ml , 25 ml, 50 ml, 100 ml, 250 ml, 500 ml, 1000 ml, 2000 ml