

Measuring Jugs



Product Categories: [Laboratory Plastic Ware](#), [Liquid-Measurement](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/measuring-jugs/>

Product Description

Measuring Jugs - Scientific Laboratory Equipment

The Measuring Jugs are moulded in Polypropylene, are clear, autoclavable and have good chemical resistance. The Measuring Jugs are provided with a handle for easy & convenient lifting and have raised graduations for easy reading. Adding to the existing range of Polylab Jugs, these short form transparent jugs, made of polypropylene are autoclavable and have non-drip spout. These jugs have raised graduations, provide excellent contact clarity, & are provided with excellent thumb grip on handle.

Features:

Scratch resistance

Sturdiness

Light weight

Particulars	Least Count	Packing
250 ml	10 ml	12 Pcs
500 ml	10 ml	12 Pcs

1000 ml	10 ml	6 Pcs
2000 ml	20 ml	6 Pcs