


# Mechanical Stage



**Product Categories:** [Microscope](#), [Microscope Assecories](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/mechanical-stage/>

## Product Description

Mechanical Stage - Laboratory Equipment

Both horizontal and vertical movements are made by a rack and pinion arrangement provided. It has a movement of 70 mm in X-axis and 30 mm in Y-axis. Its vernier reads to 0.1 mm.

**Mechanical Stage Coaxial :**

In this model both X & Y movements can be made from knobs provided on the same X-axis. The horizontal movement is made by a lead screw guide nut arrangement and vertical movement is made by a rack and pinion arrangement. It has a movement of 70 mm in X-axis and 30 mm in Y-axis. Its vernier reads to 0.1 mm