

Mecker Burner



Product Categories: [Chemistry Equipments](#), [Lab Supplies](#)

Product Tags: [chemistry lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/mecker-burner/>

Product Description

Mecker Burner - Laboratory Equipment

A Mecker burner is a laboratory burner that produces multiple open gas flames, used for heating, sterilization, and combustion. The flame burns without noise. For use with liquid Propane Gas. Base with tapered refilled inlet tube of OD 9mm.

Chrome-plated burner tube has perforated head of 25mm. Rotatable air-regulator around the jet. Overall height 16cm.

Available in size - 3/4", 1"

Feature -

Brass lacquered.

With ceramic grid cap.

With stop cock

Provide and evenly accurate flame.