

## Medical Microscope



**Product Categories:** [Microscope](#), [Microscopes](#)

**Product Tags:** [biology lab equipment](#)

**Product Page:** <https://www.labappara.com/product/medical-microscope/>

### Product Description

Medical Microscope - Laboratory Equipment

**Stand :** Stable & Robust, traditional Horse Shoe base body. Body inclinable through 90°, Standard tube length of 160 mm

**Focusing :** Separate coarse and fine focusing knobs with graduation on one side. Coarse movement by Rack & Pinion and Fine movement by Screw Lever Mechanism

**Nosepiece :** Quadruple Revolving Nosepiece is provided with positive click stop for proper centering

**Stage :** Fixed Stage size 125 X 125 mm with detachable mechanical stage for X and movement of slide up to 25 X 75 mm with graduated scales and verniers

**Illumination :** illumination through a 50mm diameter planoconcave mirror mounted on gimbals

**Objectives :** Three Achromatic Objectives : 10 X 40 X (SL) & 100 X (SL) Oil

immersion

**Eyepieces** : Two Huygenian Eyepieces : 5 X & 10 X

**Condensor** : Sub Stage Abbe condenser with N.A. 1.25 moves up and down by rack & pinion. Irs diapharagm and filter holder are mounted with the condenser

**Safety Device** : Adjustable down movement stopper for preventing accidental damage