

Melting Point Apparatus Microprocessor Controlled



Product Categories: [Lab Measuring Instruments](#)

Product Tags: [Lab Measuring Instrument](#), [Laboratory Equipment](#)

Product Page:

<https://www.labappara.com/product/melting-point-apparatus-microprocessor-controlled/>

Product Description

Melting Point Apparatus Microprocessor Controlled - Laboratory Equipment

4 digit Display. Visualization

Technique Suitable for Dark Colored/any colored samples.

Silicon Oil Bath Heating Type, Temp. Range: Up to 300°C, Suitable for Melting/Boiling Point, Readability in °C, Resolution 1°C, Programmable, Visual provision by Optical Magnifier, Automatic Ramp Rate Proceeding, Auto Gradient Rate of Heating Control, Sample Capacity.