

Melting Point Apparatus Microprocessor Controlled



Product Categories: [Analytical Lab Instruments](#)

Product Tags: [Analytical Laboratory Equipment](#)

Product Page:

<https://www.labappara.com/product/melting-point-apparatus-microprocessorcontrolled/>

Product Description

Melting Point Apparatus Microprocessor Controlled - Analytical Laboratory Equipment

LED Four digit Display.

Visualization Technique Suitable for Dark Colored/any colored samples.

Silicon Oil Bath Heating Type, Temp. Range: Up to 300°C, Suitable for Melting/Boiling Point

Readability in °C, Resolution 1°C, Programmable, Visual provision by Optical Magnifier, Automatic Ramp Rate Proceeding, Auto Gradient Rate of Heating Control, Sample Capacity.

a) Glass Beaker 200 ml, with oil level marking.

b) Oil immersion Heater Assembly complete with Teflon lid and connecting cable.

c) Magnetic capsule (35 x 9 mm) .

d) Capillary set (set of 100 nos.) Superior Quality one end Sealed.

e) Silicon Oil (200 ml) .

f) Boiling Point Tubes (Set of 5).