

# Microprocessor Flame Photometer



**Product Categories:** [Lab Measuring Instruments](#)

**Product Tags:** [Lab Measuring Equipment](#), [Laboratory Equipments](#)

**Product Page:** <https://www.labappara.com/product/microprocessor-flame-photometer/>

## Product Description

Microprocessor Flame Photometer - Laboratory Equipment

With Na & K Filter Upto 5 point calibration direct readout in ppm & meq/L, auto filter selection, auto gas cut auto ignition, large LCD Display 24x4 lines graphical through RS232 interface, single inspiration measurement results stored in memory 0 to 100% Sensitivity for full scale. Compressor & other accessories. Ca & Li filters can be provided on request at extra cost.